

MARKET-SPECIFIC INFORMATION

COSMETIC

	Anaerobes	Environmental	General Purpose Culture Media	Gram (-) Enteric Bacteria	Lactobacillus	Pseudomonas	Staphylococcus	Sterility Testing	Streptococcus	Yeast & Mold
Baird Parker Agar (7112)							X			
Bile Esculin Agar (7249)								X		
Brain-Heart Infusion Agar (7115)			X							X
Brain-Heart Infusion Broth (7116)	X		X							X
Cetrimide Agar Base (7222)						X				
Cetrimide Agar (Agar Medium N) (7688)						X				
Cooked Meat Medium (7110)	X									
D/E Neutralizing Agar (7375)		X								
D/E Neutralizing Broth (7562)		X								
D/E Neutralizing Broth with Tween (7705)		X								
Eosin Methylene Blue Agar (Holt, Harris & Teague) (7134)				X						
Eosin Methylene Blue Agar, Levine (7103)				X						
Fluid Thioglycollate Medium (7137)	X		X					X		
HC Agar Base (7520)										X
Lactobacillus MRS Agar (7543)					X					
Lactobacillus MRS Broth (7406)					X					
Lactobacillus Selective Agar Base (7234)					X					
Letheen Agar Base (7118)		X								
Letheen Agar Base, Modified (7495)		X								
Letheen Agar with Tween (7710)		X								
Letheen Agar with Tween, Modified (7712)		X								
Letheen Broth Base (7105)		X								
Letheen Broth Base, Modified (7496)		X								
Littman Agar (7173)										X
MacConkey Agar (7102)				X						
Malt Agar (7456)										X
Mannitol Salt Agar (7143)							X			
Mycobiotic Agar (7419)										X
Mycological Agar (7309)										X
Phenylethanol Agar (7147)								X		
Potato Dextrose Agar (7149)										X
Potato Dextrose Broth (7585)										X
Pseudomonas Cepacia Agar (7458)						X				
Pseudomonas Isolation Agar (7329)						X				
Rose Bengal Chloramphenicol Agar (7664)										X
Sabouraud Dextrose Agar (7150)										X
Sabouraud Dextrose Agar with Chloramphenicol (7306)										X
Sabouraud Dextrose Agar with Lecithin & Tween (7392)		X								X
Sabouraud Dextrose Broth (7617)										X
Standard Methods Agar (7157)			X							
TAT Broth (7219)								X		
Thioglycollate Medium without Indicator (7160)	X							X		
Todd Hewitt Broth (7161)	X									
Tomato Juice Agar (7349)					X					
Tryptic Soy Agar (7100)			X							
Tryptic Soy Agar with Lecithin & Tween 80 (7163)		X								
Tryptic Soy Broth (7164)			X					X		
Vogel & Johnson Agar (7207)							X			